



Metrosep A Supp 3 - 250/4,6

6.1005.320

La colonne Metrosep A Supp 3 - 250/4,6 permet de solutionner des problèmes de séparation en milieux aqueux et organiques. Elle peut être utilisée avec de nombreux éluants de manière sûre, même si la teneur en solvants organiques est importante. La colonne Metrosep A Supp 3 - 250/4.6 permet d'analyser, de manière routinière, des échantillons représentant un défi analytique, tels que la mesure d'échantillons biologiques ou la détermination des anions inorganiques dans des matrices organiques. Un gradient d'hydroxyde de sodium, donne la certitude de séparer les polyphosphates sur la colonne Metrosep A Supp 3 - 250/4,6. En mode isocratique, elle est également très bien adaptée à la séparation des sulfites, sulfates et thiosulfates en moins de 20 minutes.

Accessoires optionnels

Order no.	Description
-----------	-------------

6.1011.030

Metrosep RP 2 Guard/3.5

The Metrosep RP 2 Guard/3.5 is a guard column for universal use. It reliably protects the analytical separation column against contaminations. It reliably protects the analytical separation column against contamination, removing the smallest particles, traces of iron oxide and bacteria. The Metrosep RP 2 Guard/3.5 helps to reduce costs; its filter disk can be replaced very easily.



6.1011.040

Metrosep RP 3 Guard HC/4.0

The Metrosep RP 3 Guard HC/4.0 is a guard column for universal use. It reliably protects the analytical separation column against contamination, securely removing lipophilic organic contamination, the smallest particles, traces of iron oxide and bacteria. The guard column is based on a polymer material and thanks to its larger pack volume, has a significantly higher capacity than the Metrosep RP 2 Guard/3.5. It can be used throughout the entire pH range.



6.1011.130

Spare filter for RP 2 Guard/3.5, 10 pieces

