



Metrosep A Supp 5 - 250/2,0

6.1006.230

La Metrosep A Supp 5 - 250/2,0 est la colonne de séparation microbore haute performance, qui permet de résoudre facilement et de manière reproductible même les problèmes de séparation complexes. Le spectre des applications possibles de cette colonne dépasse largement la seule détection des anions standard. La Metrosep A Supp 5 - 250/2,0 est utilisée chaque fois qu'il s'agit d'associer des performances de séparation extrêmes à des limites de sensibilité très basses et une très faible consommation d'éluant.

Cette colonne convient particulièrement au couplage CI-SM du fait du faible débit d'éluant.

Accessoires optionnels

Order no.	Description
-----------	-------------

The Metrosep A Supp 5 Guard/2.0 reliably protects the 2 mm versions of the Metrosep A Supp 5 and Metrosep A Supp 7 IC anion columns against contamination from the sample or eluent.

It contains the same separation material as the Metrosep A Supp 5, is also made of PEEK, and is screwed directly onto the respective separation column with virtually no dead volume ("On Column Guard System"). The guard column prolongs the service life of the analytical column, with practically no influence on its chromatographic separating efficiency. The economical price and simple handling make using the A Supp 5 Guard/2.0 highly recommended.



The Metrosep A Supp 5 S-Guard/2.0 reliably protects the 2 mm versions of the Metrosep A Supp 5 and Metrosep A Supp 7 IC anion columns against contamination from the sample or eluent.

It contains the same separation material as the Metrosep A Supp 5, is also made of PEEK, and is connected directly onto the respective separation column via capillary connection. The guard column prolongs the service life of the analytical column, with practically no influence on its chromatographic separating efficiency. The economical price and simple handling make using the A Supp 5 S-Guard/2.0 highly recommended.

