



# Metrosep A Supp 5 - 100/4,0

6.1006.510

La Metrosep A Supp 5 - 100/4,0 permet des séparations hautement efficaces et rapides. Cette propriété a permis à la Metrosep A Supp 5 - 100/4,0 de s'établir comme colonne standard pour des analyses rapides et la détermination des anions à élution tardive comme le perchlorate.

## Accessoires optionnels

| Order no.  | Description                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 6.1006.500 | Metrosep A Supp 5 Guard/4.0                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|            | <p>The Metrosep A Supp 5 Guard/4.0 reliably protects the Metrosep A Supp 5 and 7 IC anion columns against contamination from the sample or the eluent.</p> <p>It contains the same separation material as the Metrosep A Supp 5, is also made of PEEK, and is screwed directly onto the respective separation column with virtually no dead volume ("On Column Guard System"). The guard column prolongs the service life of the analytical column, with practically no influence on its chromatographic separating efficiency. The economical price and simple handling make using the A Supp 5 Guard/4.0 highly recommended.</p> |



6.1006.540

#### Metrosep A Supp 5 S-Guard/4.0

The Metrosep A Supp 5 S-Guard/4.0 reliably protects the Metrosep A Supp 5 and 7 IC anion columns against contaminations from the sample or eluent.



---

6.1020.520

#### Metrosep A Supp 10 Guard HC/4.0

The Metrosep A Supp 10 Guard HC/4.0 is the high-capacity version of the Metrosep A Supp 10 Guard/4.0.

