



Metrosep A Supp 5 - 50/4,0

6.1006.550

La Metrosep A Supp 5 - 50/4,0 sépare les sept anions standard en moins de 6 minutes. Même le fluorure reste séparé du pic de temps mort et se laisse parfaitement intégrer. Cette colonne basée sur un polymère à base d'alcool polyvinylique se distingue, comme toutes les autres colonnes de la famille A Supp 5, par un nombre élevé de plateaux et, par conséquent, par une excellente performance de séparation. La Metrosep A Supp 5 - 50/4,0 est la colonne de choix lorsqu'il s'agit de résoudre de simples tâches de séparation en peu de temps – sans renoncer à des limites de sensibilité très basses.

Accessoires optionnels

Order no.	Description
-----------	-------------

The Metrosep A Supp 5 Guard/4.0 reliably protects the Metrosep A Supp 5 and 7 IC anion columns against contamination from the sample or the eluent.

It contains the same separation material as the Metrosep A Supp 5, is also made of PEEK, and is screwed directly onto the respective separation column with virtually no dead volume ("On Column Guard System"). The guard column prolongs the service life of the analytical column, with practically no influence on its chromatographic separating efficiency. The economical price and simple handling make using the A Supp 5 Guard/4.0 highly recommended.



The Metrosep A Supp 5 S-Guard/4.0 reliably protects the Metrosep A Supp 5 and 7 IC anion columns against contaminations from the sample or eluent.

