



Metrosep C 4 - 100/4.0

6.1050.410

La version 100 mm de la colonne Metrosep C 4 est conçue pour les déterminations rapides des cations standard. On obtient des temps de rétention très courts, les temps d'élution du sodium et de l'ammonium se différenciant tout de même de 25 secondes. En utilisant un éluant spécial, la colonne Metrosep C 4 - 100/4,0 permet de déterminer en moins de 5 minutes six cations : lithium, ammonium, sodium, calcium, magnésium et potassium.

Contenu de la livraison 6.1050.410

Qt.	Order no.	Description
2 PCS	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Accessoires optionnels

Order no.	Description	
6.1050.500	Metrosep C 4 Guard/4.0 The Metrosep C 4 Guard/4.0 contains the C-4 column material and is used to protect all Metrosep cation columns that have a substrate based on silica gel. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column. The economical price is an additional plus. The Metrosep C 4 Guard/4.0 also works based on the "On Column Guard System" and is screwed directly onto the separation column with virtually no dead volume.	 A photograph of the Metrosep C 4 Guard/4.0 column. It is a small, cylindrical, tan-colored component with white plastic end caps. A small label is attached to the side, showing the product name and a diagram of the column's internal structure.
6.1050.510	Metrosep C 4 S-Guard/4.0 The Metrosep C 4 S-Guard/4.0 contains the C-4 column material and is used to protect all Metrosep cation columns that have a substrate based on silica gel. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column. The economical price is an additional plus.	 A photograph of the Metrosep C 4 S-Guard/4.0 column. It is a small, cylindrical, tan-colored component with white plastic end caps. A small label is attached to the side, showing the product name and a diagram of the column's internal structure.
6.1050.530	Metrosep C 4 S-Guard - 50/4.0 The Metrosep C 4 S-Guard - 50/4.0 contains the C 4 column material and is used to protect all Metrosep cation columns that have a substrate based on silica gel. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column. In comparison with Standard C 4 guard columns, the Metrosep C 4 S-Guard - 50/4.0 has a greater capacity and therefore an even longer service life.	 A photograph of the Metrosep C 4 S-Guard - 50/4.0 column. It is a small, cylindrical, tan-colored component with green plastic end caps. A small label is attached to the side, showing the product name and a diagram of the column's internal structure.