



Metrosep C 4 - 250/4.0

6.1050.430

La Metrosep C 4 - 250/4,0 est la colonne pour cations avec la plus haute capacité de la série C 4. Elle est destinée aux applications exigeant des performances de séparation très élevées. Cette colonne permet d'analyser de manière fiable des échantillons présentant des différences de concentration extrêmes. Les performances de la colonne se dévoilent par exemple lors de l'analyse de l'eau d'alimentation de chaudières, où il s'agit de quantifier correctement 7 µg/L de sodium en présence de 7 mg/L de monoéthanolamine (MEA). Avec la C 4 - 250/4,0, à la fois les amines et les métaux de transition ainsi que les métaux alcalins et alcalino-terreux peuvent être déterminés en une seule opération.

Contenu de la livraison 6.1050.430

Qt.	Order no.	Description
-----	-----------	-------------

2 PCS

6.2744.060

Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Accessoires optionnels

Order no.	Description	
6.1050.500	Metrosep C 4 Guard/4.0	
	The Metrosep C 4 Guard/4.0 contains the C-4 column material and is used to protect all Metrosep cation columns that have a substrate based on silica gel. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column. The economical price is an additional plus. The Metrosep C 4 Guard/4.0 also works based on the "On Column Guard System" and is screwed directly onto the separation column with virtually no dead volume.	
6.1050.510	Metrosep C 4 S-Guard/4.0	
	The Metrosep C 4 S-Guard/4.0 contains the C-4 column material and is used to protect all Metrosep cation columns that have a substrate based on silica gel. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column. The economical price is an additional plus.	
6.1050.530	Metrosep C 4 S-Guard - 50/4.0	
	The Metrosep C 4 S-Guard - 50/4.0 contains the C 4 column material and is used to protect all Metrosep cation columns that have a substrate based on silica gel. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column. In comparison with Standard C 4 guard columns, the Metrosep C 4 S-Guard - 50/4.0 has a greater capacity and therefore an even longer service life.	