



Rack d'échantillons, 48 x 75 mL (PVC)

6.2041.350

Rack d'échantillons pour 815 Robotic USB Sample Processor XL ou 855 Robotic Titrosampler « Basic » pour 48 béchers d'échantillon. Vous pouvez utiliser les béchers d'échantillon suivants : 6.1432.210, 6.1459.400.

Plastique : polychlorure de vinyle (PVC)

Accessoires optionnels

Order no.	Description
6.1432.210	Sample beaker, clear glass, 75 mL, 1 piece
	Sample beaker, clear glass, 75 mL, suitable for sample racks: 6.2041.340, 6.2041.350, 6.2041.800 and 6.2041.900, 1 piece



6.1458.070

Titration head / 2 x SGJ 14

Titration head for 75 mL beakers. With 2 openings SGJ 14, 2 openings 6.4 mm, and 3 pre-mounted dosing tips. Suitable for rinsing and aspiration equipment.



6.1459.400

Sample beaker, plastic (PE), 75 mL, 144 pieces

Sample beaker, plastic (polyethylene, PE) 75 mL, suitable for sample racks: 6.2041.340, 6.2041.350, and 6.2041.800, suitable for 6.02041.040 OMNIS Software sample rack, 144 pieces.



6.1462.050

Robotic arm with titration head for 786 Swing Head, left or right swinging

Titration head for Swing Head

