



Crépine d'aspiration

6.2821.090

Dimensions des pores 27 µm. Jeu de 5 pièces. Pour tuyau d'aspiration (6.1834.000) et tubes de filtration (6.1821.040) et (6.1821.050).

Accessoires optionnels

Order no.	Description
6.02744.050	Adapter Luer inner / M6 outer

Adapter for syringe for pre-rinsing the aspiration filter (6.2821.090).
The adapter can be screwed onto the aspiration filter. A syringe can be used to rinse the aspiration filter with ultrapure water or eluent.



6.2816.020

Syringe 10 mL with Luer connection

For various applications in IC and VA

