



Metrosep Amino Acids 1 - 100/4,0

6.4001.410

Colonne de séparation pour la détermination des acides aminés au moyen de la détection UV/VIS après PCR à la ninhydrine à 120 °C. La Metrosep Amino Acids 1 - 100/4,0 est la colonne de séparation standard des aminoacides. Cette colonne est basée sur un matériau de polystyrène-divinylbenzène sulfoné. La détermination des acides aminés s'effectue au moyen de la détection photométrique après dérivatisation post-colonne avec de la ninhydrine.

La Metrosep Amino Acids 1 - 100/4.0 permet une séparation d'un maximum de 44 acides aminés dont tous les acides aminés survenant naturellement dans la recherche et la pratique.

Contenu de la livraison 6.4001.410

Qt.	Order no.	Description
2 PCS	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Accessoires optionnels

Order no.	Description	
6.1011.030	Metrosep RP 2 Guard/3.5	
	The Metrosep RP 2 Guard/3.5 is a guard column for universal use. It reliably protects the analytical separation column against contaminations. It reliably protects the analytical separation column against contamination, removing the smallest particles, traces of iron oxide and bacteria. The Metrosep RP 2 Guard/3.5 helps to reduce costs; its filter disk can be replaced very easily.	
6.1011.040	Metrosep RP 3 Guard HC/4.0	
	The Metrosep RP 3 Guard HC/4.0 is a guard column for universal use. It reliably protects the analytical separation column against contamination, securely removing lipophilic organic contamination, the smallest particles, traces of iron oxide and bacteria. The guard column is based on a polymer material and thanks to its larger pack volume, has a significantly higher capacity than the Metrosep RP 2Guard/3.5. It can be used throughout the entire pH range.	
6.1011.130	Spare filter for RP 2 Guard/3.5, 10 pieces	