



Équipement pour la détermination de la correction de température pour les Rancimat et les thermostats PVC.

6.5616.100

Set permettant un ajustement exact de la température

Contenu de la livraison 6.5616.100

Qt.	Order no.	Description
1 PCS	6.1111.010	Pt100 temperature sensor for stability measuring devices
Pt100 temperature sensor (length 17.5 cm) with fixed cable and mini-DIN-plug. This glass-free temperature sensor is well suited for use with a 711 Liquino and instruments for stability measurement (743 Rancimat, 763 PVC Thermomat, 892 Professional Rancimat, or 895 Professional PVC Thermomat).		



2 PCS

6.2042.040

Clamp for temperature sensor

Holding clamp for Pt100 temperature sensor for stability measuring instruments



1 PCS

6.2326.000

Silicone oil for stability measuring instruments (50 mL)

Silicone oil for the determination of temperature correction in stability measuring instruments



1 PCS

6.2418.000

Air tube for oil measurements

Set of 12 pieces. In combination with 6.1429.040.

