



Kit de recharge pour Lumitrode O2

6.5623.000

Kit de recharge pour Lumitrode O₂. Le kit comprend un capuchon O₂ homologué de recharge, un récipient de calibrage et 5 étalons oxygène de 30 mL à 0 %.

Contenu de la livraison 6.5623.000

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS	6.1260.000	Calibration vessel O2
-------	------------	-----------------------

Calibration vessel for calibrating the O₂-Lumitrode at 100% air saturation.



Oxygen standard for 0% calibration of the O₂-Lumitrode. Contains 5 sachets of 30 mL each.



Accessoires optionnels

Order no.	Description
6.1116.000	O2 Lumitrode The optical sensor for measuring dissolved oxygen (DO) can be used with a 913 pH/DO meter or with a 914 pH/DO conductometer. The measuring principle of the sensor is based on luminescence quenching. The space-saving and maintenance-free sensor is suitable for DO measurement, for example in: • Water quality control • Wastewater industry • Beverage production • Fish farming



This sensor is supplied with a calibration vessel.

Where necessary, it is easy to replace the measurement cap (O₂ cap) that contains the oxygen-intensive luminophore.
