



TacticID-1064 ST

BWT-840001134

Le TacticID®-1064 ST est un appareil d'analyse portable Raman 1 064 nm pour une identification rapide sur le terrain d'explosifs, de stupéfiants et d'autres matières suspecte. La capacité « See-through » (transparence) du TacticID-1064 ST peut analyser des échantillons de manière non destructive à travers un emballage opaque et transparent, avec un niveau de menace d'échantillon affiché en évidence pour les premiers intervenants, le personnel de sécurité, les forces de l'ordre, les démineurs, les douanes et la protection des frontières, et les équipes de gestion de matières dangereuses pour leur permettre d'agir rapidement avec un minimum de contact avec l'échantillon.

Le TacticID-1064 ST utilise la spectroscopie Raman éprouvée, en combinaison avec la technologie brevetée STRaman®, permettant aux utilisateurs d'obtenir une identification exploitable en temps réel de produits chimiques inconnus, de stupéfiants, de médicaments, d'explosifs et de nombreuses autres substances, même à travers des barrières opaques, réduisant considérablement l'incertitude opérationnelle et le temps de réponse.

Le TacticID-1064 ST, avec une excitation laser à 1 064 nm et un adaptateur ST pour les applications « see-through » (transparence), exploite une vaste gamme d'échantillons, délivrant un spectre sans fluorescence, sans brûlure d'échantillon, permettant aux utilisateurs d'identifier les matériaux sombres et colorés, les échantillons de rue difficiles, les mélanges et les matériaux hétérogènes directement à travers l'emballage.

Ce système IP68 dispose d'un écran haute luminosité avec écran tactile et/ou interface matérielle par boutons pour faciliter l'utilisation, même à travers un équipement de protection.

Contenu de la livraison BWT-840001134

Qt.	Order no.	Description	
1 PCS	BWT-840000102	Laser safety goggles (1064 nm) Laser safety goggles for lasers with excitation at 1,064 nm – suitable for Class IIIB laser systems from B&W Tek.	
1 PCS	BWT-840000295	Vials for liquids, 15 mm in diameter Borosilicate glass vials (diameter 15 mm), package of 6, with pouch.	
1 PCS	BWT-840000329	Pouch for 15 mm Vials Pouch for 15 mm vials.	
1 PCS	BWT-840000344	Vial Holder Adapter Vial holder adapter, compatible with 15 mm-diameter vials.	

1 PCS

BWT-840000345 Handheld Raman bottle adapter

Handheld Raman accessories: bottle adapter

NR2-BSA



1 PCS

BWT-840000348 Polystyrene Validation Cap


Polystyrene Validation cap with integrated ASTM standard polystyrene.



1 PCS

BWT-840000353 Shaft Protection Cap

Handheld Raman shaft protective cap, package of 10.



1 PCS

BWT-840000419 TacticID Stylus/pen combo

Stylus/pen combo for touch screens (package of 5).

TID-STY



1 PCS

BWT-840000858 USB cable

USB cable supporting micro USB port up to type A USB 2.0.



1 PCS

BWT-840000948 USB OTG Host Cable

USB OTG host cable for data transmission between a USB flash drive and a handheld system.



1 PCS

BWT-840000949 Large Spot Adapter

The large-spot sampling adapter is designed for use with Metrohm handheld Raman systems, providing a larger spot size and increased sampling depth to minimize the risk of sample burning, particularly for sensitive or dark-colored samples. Additionally, this adapter supports see-through capabilities when used with devices compatible with see-through software.



TID-STA

Accessoires optionnels

Order no.	Description	
BWT-840000294	Battery Charging Cradle	
	Lithium-ion battery charging cradle.(for charging NR-BAT and TID-BAT only).	
BWT-840000361	Handheld Laser Power Meter	
	Laser Power Meter, Pen Type, Calibrated. Range: 5 micro-Watts to 1 Watt. Used in OQ/PQ and annual certification laser power measurement.	
BWT-840000414	18 V Power Adapter	
	Power adapter for 18 V DC, 1.67 A. Plugs for USA, EU, UK and AU included in the scope of delivery.	
	TID-PWA	
BWT-840000428	Wireless Router	
	Wireless, battery-powered portable router; 3G/4G connectivity.	

BWT-840000498 Disposable Battery Cartridge

TacticID battery cartridge for use with six type 123A disposable batteries (e.g. Surefire). Batteries not included.



BWT-840000969 Handheld Raman 1064 Immersion Probe

Handheld Raman Spectrometer 1064 accessory:. Contact immersion probe, 6" length, working distance approx. 1 mm, operating pressure up to 2 bar at room temperature.

NRX-CTP-IMP



BWT-840000970 Handheld Raman large bottle adapter (1064)

Handheld Raman Spectrometer 1064 accessories: Large bottle adapter with 10 mm working distance.

NRX-BSA



BWT-840001154 Shoulder bag for HH Raman Series

Shoulder bag for TacticID series. Constructed with 1000d heavy duty fabric, multiple tactical pouches and compartment designed for rugged use and hard wearing.. Compatible with all TacticID products and NanoRam products. For transportation purposes only.

TID-SBG

