



Stampante termica Custom Q3X

2.142.0100

Stampante compatta con interfaccia USB per

- 900 Touch Control
- 915 KF Ti-Touch
- 916 Ti-Touch
- 917 Coulometer
- 877 / 848 Titrino plus
- 865 / 876 Dosimat plus
- 91X Meter (cavo 6.2151.140)
- Eco Dosimat / Titrator
- 862 Compact Titrosampler
- 870 KF Titrino plus
- 899 Coulometer

Larghezza carta 60 mm (40 caratteri). Cavo USB 6.2151.120 incluso.

Parti incluse 2.142.0100

Qt.	Order no.	Descrizione
-----	-----------	-------------

1 PCS	1.142.0100	Custom Q3X thermal printer
-------	------------	----------------------------

Compact printer with USB interface for

- 900 Touch Control
- 915 KF Ti-Touch
- 916 Ti-Touch
- 917 Coulometer
- 877 / 848 Titrino plus
- 865 / 876 Dosimat plus
- 91X Meter
- Eco Dosimat / Titrator



Paper width 60 mm (40 characters). Including 6.2151.120 USB cable.

1 PCS	6.2151.120	Cable USB MINI (OTG) - USB B, 1.8 m
-------	------------	-------------------------------------

For connecting USB instruments.



Accessori opzionali

Order no.	Descrizione	
6.2151.140	Y cable USB A St - USB B St - Mini B St Y cable for connection of a USB printer to the pH/Conductometers 912 / 913 / 914. This cable allows a printer and the power supply unit to be connected to the measuring instrument at the same time.	
6.9001.050	Thermal paper for Custom Neos-N / Q3X Thermal paper 60 mm / 25 m for Custom Neos-N / Q3X. New: 25-year legibility with correct storage. 1 roll	