



728 Magnetic Stirrer con stand

2.728.0040

Per l'estensione del 645, 685, 665/765, 715/775, 725/776 Dosimat al Titration Stands, o completare per Titrino e KF Coulometer; incl. piastra di base 6.2001.010 con asta di supporto e portaelettrodo; senza adattatore, con cavo di connessione fisso.

Parti incluse 2.728.0040

Qt.	Order no.	Descrizione
1 PCS	1.728.0010	728 Stirrer without stand For mounting on a stand rod 10 mm dia.; without base plate (or holder) and without stand rod; without power adapter, with fixed connecting cable.



1 PCS

6.1903.020

Stirring bar / 16 mm

Stirring Bar with magnetic core, PTFE coated, length 16 mm



1 PCS

6.1903.030

Stirring bar / 25 mm

Stirring bar with magnetic core, PTFE covering, length 25 mm.



1 PCS

6.2001.010

Support stand with support rod

For Dosimats and Titrinos



1 PCS

6.2013.010

Clamping ring

For support rods with a diameter of 10 mm.



1 PCS

6.2021.020

Electrode holder

Electrode holder for 4 electrodes and 2 buret tips



1 PCS

6.2723.130

Cover for 655/665/765

Dust cover for 655, 665, 765 Dosimats

