



827 pH lab; 230 V, EU con Primatrode

2.827.0114

pH-metro da laboratorio con interfaccia IrDa e Primatrode 6.0228.010

Parti incluse 2.827.0114

Qt.	Order no.	Descrizione
1 PCS	6.2013.010	Clamping ring

For support rods with a diameter of 10 mm.



1 PCS

6.2016.050

Support rod / 300 mm



1 PCS

6.2021.020

Electrode holder

Electrode holder for 4 electrodes and 2 buret tips



1 PCS

6.2161.010

Plug-in power supply unit EU 230V /6V DC

Power supply for 827 pH lab



1 PCS

6.2621.070

5 mm hex key for IC Sample Processors



Accessori opzionali

Order no.	Descrizione	
6.2307.230	Buffer solutions pH 4, 7 and 9 Mixed buffer solutions pH 4.00/7.00/9.00 (25 °C) in single use sachets, colourless, box of 3 x 10 x 30 mL	
6.2325.000	pHit kit Maintenance kit for pH electrodes The kit contains: • 50 mL cleaning solution • 50 mL 3M KCl solution • 50 mL storage solution • 2 Storage vessels • Instructions for use	