



849 Level Control

2.849.0010

Kit di accessori aggiuntivi per Sample Processors per controllare il livello di riempimento del recipiente di lavaggio o di scarico tramite Remote. Per evitare che le pompe lavorino a secco e/o la fuoriuscita dal contenitore dello scarico. Adatto per essere usato con soluzioni acquose, solventi e sospensioni.

Parti incluse 2.849.0010

Qt.	Order no.	Descrizione
1 PCS	1.849.0020	849 Level Control metal-free (for canisters)

Additional equipment for Metrohm instruments. For monitoring the filling level of bottles for eluent, auxiliary or waste solutions via remote control. The 849 Level Control prevents pumps from running dry and/or canisters from overflowing and is suitable for use with aqueous solutions, solvents and suspensions.

The metal-free sensor works in solutions with very low conductivity. Therefore it may be used in ion chromatography.



1 PCS

6.1113.000

Sensor for 849 Level Control

Sensor for 849 Level Control



1 PCS

6.2151.060

Cable for Level Control and Eluent Production Module

Connecting cable for sensor and Level Control or Eluent Production Module.

