



## Micro Au Titrode (lunghezza di inserimento 30,8 cm)

6.00435.120

Elettrodo ad anello d'oro combinato con una membrana di vetro per pH come elettrodo di riferimento, lunghezza di inserimento = 30,8 cm, rastremato a un diametro di 6,4 mm.

Questo elettrodo privo di manutenzione è adatto alle titolazioni redox con valore di pH costante, ad es.:

- Determinazione della domanda chimica di ossigeno (COD, ferrometria)
- Determinazione di penicillina e ampicillina
- Titolazioni con  $\text{Hg}(\text{NO}_3)_2$
- Titolazioni in presenza di cromo e ferro

Questo elettrodo viene conservato in acqua distillata.

Parti incluse 6.00435.120

Qt.	Order no.	Descrizione
-----	-----------	-------------

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



---

1 PCS


6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



# Accessori opzionali

Order no.	Descrizione	
6.2104.020	Electrode cable / 1 m / F  For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).	

---