



## iAg Titrode con rivestimento AgBr

6.00470.300BR

Elettrodo ad anello d'argento combinato, intelligente con membrana di vetro per pH come elettrodo di riferimento e chip di memorizzazione integrato per i dati del sensore.

L'anello d'argento è rivestito con bromuro (AgBr) a garanzia di una maggiore sensibilità e migliori limiti di prova.

Questo elettrodo, privo di manutenzione, è adatto alle titolazioni per precipitazione con valore di pH costante (titolante nitrato di argento) ad es. di:

- cloruro, bromuro, ioduro
- solfuri
- idrogeno solforato
- mercaptano
- cianuro

Questo elettrodo viene conservato in acqua distillata.

Gli elettrodi iTrode possono essere utilizzati con Titrando, Ti-Touch o con i misuratori 913/914.

Parti incluse 6.00470.300BR

| Qt.                                                                                 | Order no.  | Descrizione                                                                                                                    |
|-------------------------------------------------------------------------------------|------------|--------------------------------------------------------------------------------------------------------------------------------|
| 1 PCS                                                                               | 6.1236.050 | Sleeve with SGJ 14/12 mm<br>Sleeve with SGJ 14/12 mm, polyethylene.                                                            |
|  |            |                                                                                                                                |
| 1 PCS                                                                               | 6.2008.040 | Storage vessel<br>Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units. |
|  |            |                                                                                                                                |

# Accessori opzionali

| Order no.  | Descrizione                                                                                                                                                                                                                                              |
|------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 2.854.0010 | <div>854 iConnect</div> <div>The 854 iConnect – Electrode cable and measuring amplifier for intelligent "iTrodes" electrodes (cable length 150 cm).</div> <div></div> |

---