



## iAg Titrode con rivestimento Ag<sub>2</sub>S

6.00470.300S

Elettrodo ad anello d'argento combinato, intelligente con membrana di vetro per pH come elettrodo di riferimento e chip di memorizzazione integrato per i dati del sensore.

L'anello d'argento è rivestito con solfuro (Ag<sub>2</sub>S) a garanzia di una maggiore sensibilità e migliori limiti di prova.



Questo elettrodo, privo di manutenzione, è adatto alle titolazioni per precipitazione con valore di pH costante (titolante nitrato di argento) ad es. di:

- cloruro, bromuro, ioduro
- solfuri
- idrogeno solforato
- mercaptano
- cianuro

Questo elettrodo viene conservato in acqua distillata.


Gli elettrodi iTrode possono essere utilizzati con Titrando, Ti-Touch o con i misuratori 913/914.

Parti incluse 6.00470.300S

Qt.	Order no.	Descrizione
1 PCS	6.1236.050	Sleeve with SGJ 14/12 mm Sleeve with SGJ 14/12 mm, polyethylene.
		
1 PCS	6.2008.040	Storage vessel Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.
		

# Accessori opzionali

Order no.	Descrizione
2.854.0010	854 iConnect  The 854 iConnect – Electrode cable and measuring amplifier for intelligent "iTrodes" electrodes (cable length 150 cm).



---