



Metrosep A Supp 18 - 250/4,0

6.01033.430

La selettività delle colonne Metrosep A Supp 18 è specificamente progettata per lavorare con eluenti di idrossido. La Metrosep A Supp 18 - 250/4,0 è la versione più lunga tra le colonne A Supp 18. Grazie alle dimensioni ridotte delle particelle di soli 3,5 µm, si raggiungono elevate prestazioni di separazione. I picchi simmetrici consentono la facile determinazione di anioni, anche in bassi intervalli di concentrazione µg/L. La Metrosep A Supp 18 - 250/4,0 è particolarmente adatta per la separazione degli ossalogenuri clorite, bromato e clorato oltre agli anioni standard a temperatura ambiente.

Accessori opzionali

Order no.	Descrizione
-----------	-------------

The Metrosep A Supp 18 Guard/4.0 reliably protects the Metrosep A Supp 18 analytical separation column against contamination from sample or eluent. It contains the same separation material as the Metrosep A Supp 18 IC Columns, is also made of PEEK, and is screwed directly onto the respective separating column with virtually no dead volume ("On Column Guard System"). The guard column prolongs the service life of the analytical column, without influencing chromatographic separating efficiency.

