



Metrosep C Supp 2 - 100/4,0

6.01053.410

Il materiale di separazione della Metrosep C-Supp-2 è basato su un copolimero di polistirene-divinilbenzene con gruppi carbossilici. È adatta per la separazione di cationi monovalenti e divalenti. La Metrosep C Supp 2 - 100/4,0 è la colonna di separazione più corta della famiglia Metrosep C-Supp-2. È adatta in particolare all'analisi delle tracce di cationi standard. I limiti di quantificazione al di sotto dell'intervallo $\mu\text{g/L}$ vengono raggiunti grazie ai rumori della linea di riferimento estremamente basso dopo soppressione sequenziale.

Parti incluse 6.01053.410

Qt.	Order no.	Descrizione
2 PCS	6.2744.060	Threaded stopper
For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.		



Accessori opzionali

Order no.	Descrizione
6.01053.500	Metrosep C Supp 2 Guard/4.0

The Metrosep C Supp 2 Guard/4.0 contains the C Supp 2 column material and is used to protect Metrosep C Supp 2 cation columns. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column. The Metrosep C Supp 2 Guard/4.0 also works based on the "On Column Guard System" and is screwed directly onto the separation column with virtually no dead volume.

