



Metrosep A Supp 5 - 250/2,0

6.1006.230

La Metrosep A Supp 5 - 250/2,0 è la colonna di separazione microbore ad alta prestazione con cui possono essere risolti in modo facile e riproducibile anche complessi problemi di separazione. Il campo di applicazione di questa colonna supera di gran lunga l'individuazione degli anioni standard. La Metrosep A Supp 5 - 250/2,0 viene sempre utilizzata laddove occorre combinare la massima efficienza di separazione ai limiti di rilevazione più bassi e al minor consumo di eluente.

Grazie alla ridotta portata di eluente, questa colonna è particolarmente adatta per l'accoppiamento IC-MS.

Accessori opzionali

Order no.	Descrizione
-----------	-------------

The Metrosep A Supp 5 Guard/2.0 reliably protects the 2 mm versions of the Metrosep A Supp 5 and Metrosep A Supp 7 IC anion columns against contamination from the sample or eluent.

It contains the same separation material as the Metrosep A Supp 5, is also made of PEEK, and is screwed directly onto the respective separation column with virtually no dead volume ("On Column Guard System"). The guard column prolongs the service life of the analytical column, with practically no influence on its chromatographic separating efficiency. The economical price and simple handling make using the A Supp 5 Guard/2.0 highly recommended.



The Metrosep A Supp 5 S-Guard/2.0 reliably protects the 2 mm versions of the Metrosep A Supp 5 and Metrosep A Supp 7 IC anion columns against contamination from the sample or eluent.

It contains the same separation material as the Metrosep A Supp 5, is also made of PEEK, and is connected directly onto the respective separation column via capillary connection. The guard column prolongs the service life of the analytical column, with practically no influence on its chromatographic separating efficiency. The economical price and simple handling make using the A Supp 5 S-Guard/2.0 highly recommended.

