



Metrosep A Supp 5 - 50/4,0

6.1006.550

I sette anioni standard possono essere separati dalla Metrosep A Supp 5 - 50/4,0 in meno di 6 minuti. Anche il fluoruro viene separato dal picco di iniezione e può essere perfettamente integrato. La colonna che si basa su di un polimero alcool polivinilico si contraddistingue, come tutte le colonne della famiglia A Supp 5, per un numero elevato di piatti e quindi ottime prestazioni di separazione. La Metrosep A Supp 5 - 50/4,0 è la colonna di scelta quando devono essere eseguiti semplici compiti di separazione in breve tempo, senza rinunciare a limiti di determinazione molto bassi.

Accessori opzionali

Order no.	Descrizione
-----------	-------------

The Metrosep A Supp 5 Guard/4.0 reliably protects the Metrosep A Supp 5 and 7 IC anion columns against contamination from the sample or the eluent.

It contains the same separation material as the Metrosep A Supp 5, is also made of PEEK, and is screwed directly onto the respective separation column with virtually no dead volume ("On Column Guard System"). The guard column prolongs the service life of the analytical column, with practically no influence on its chromatographic separating efficiency. The economical price and simple handling make using the A Supp 5 Guard/4.0 highly recommended.



The Metrosep A Supp 5 S-Guard/4.0 reliably protects the Metrosep A Supp 5 and 7 IC anion columns against contaminations from the sample or eluent.

