



## Metrosep A Supp 16 - 150/4,0

6.1031.420

La Metrosep A Supp 16 è particolarmente adatta per problemi di separazione ad alta capacità e si contraddistingue per l'eccellente risoluzione. La Metrosep A Supp 16 - 150/4,0 è basata su un copolimero di polistirene/divinilbenzene a superficie funzionalizzata. I gruppi funzionali sono collegati in modo covalente.

La Metrosep A Supp 16 - 150/4,0 possiede una risoluzione eccellente e risolve i più complessi problemi di separazione. La colonna è particolarmente adatta per applicazioni che presentano un elevato carico ionico ma che non richiedono la massima risoluzione. È una delle colonne standard nella cromatografia anionica.

Parti incluse 6.1031.420

Qt.	Order no.	Descrizione
-----	-----------	-------------

2 PCS



6.2744.060

Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



## Accessori opzionali

Order no.	Descrizione	
6.1031.500	Metrosep A Supp 16 Guard/4.0	
	<p>The Metrosep A Supp 16 Guard/4.0 reliably protects the Metrosep A Supp 16 analytical separation columns against contamination. Thanks to the "On Column Guard System", the guard column is very easy to handle. The guard column screws easily onto the analytical column. No tools are required.</p>	
6.1031.510	Metrosep A Supp 16 S-Guard/4.0	
	<p>The Metrosep A Supp 16 S-Guard/4.0 reliably protects the Metrosep A Supp 16 analytical separation columns against contamination. The guard column is connected to the separation column simply using capillary connections. No tools are required.</p>	