



Metrosep C 4 - 150/4,0

6.1050.420

La Metrosep C 4 - 150/4,0 è la colonna standard universale per l'analisi dei cationi con elevata prestazione di separazione in breve tempo. La Metrosep C 4 - 150/4,0 è la colonna di separazione ideale per l'analisi dei metalli alcalini e alcalino-terrosi in mezzi acquosi.




Parti incluse 6.1050.420

Qt.	Order no.	Descrizione
2 PCS	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Accessori opzionali

Order no.	Descrizione	
6.1050.500	Metrosep C 4 Guard/4.0	
	<p>The Metrosep C 4 Guard/4.0 contains the C-4 column material and is used to protect all Metrosep cation columns that have a substrate based on silica gel. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column. The economical price is an additional plus. The Metrosep C 4 Guard/4.0 also works based on the "On Column Guard System" and is screwed directly onto the separation column with virtually no dead volume.</p>	
6.1050.510	Metrosep C 4 S-Guard/4.0	
	<p>The Metrosep C 4 S-Guard/4.0 contains the C-4 column material and is used to protect all Metrosep cation columns that have a substrate based on silica gel. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column. The economical price is an additional plus.</p>	
6.1050.530	Metrosep C 4 S-Guard - 50/4.0	
	<p>The Metrosep C 4 S-Guard - 50/4.0 contains the C 4 column material and is used to protect all Metrosep cation columns that have a substrate based on silica gel. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column. In comparison with Standard C 4 guard columns, the Metrosep C 4 S-Guard - 50/4.0 has a greater capacity and therefore an even longer service life.</p>	