



Tappi in PEEK GC 29 per cella di titolazione KFC


6.1446.240

Tappi per titolazione volumetrica in una cella di titolazione (coulometria)

Parti incluse 6.1446.240

Qt.	Order no.	Descrizione
1 PCS	6.01406.010	Adsorber tube with lid

Adsorber tube with lid for OMNIS Karl Fischer titrator, volumetric



1 PCS

6.02709.010

Stopper

Plastic stoppers, suitable for 2.1001.4220 OMNIS Karl Fischer Titrator.



1 PCS

6.02709.030

Screw nipple

Plastic screw nipple, suitable for 6014050X0 titration vessel lids.



1 PCS

6.1543.200

Titration tip / M6 thread

Titration tip with antidiffusion valve.



1 PCS

6.2713.010

PTFE sleeve SGJ29

Sleeve for fat-free ground-joint connection with the standard ground-joint SGJ 29.

