



## Beaker per campioni in plastica (PP) da 120 mL, 250 pezzi

6.1459.300

Beaker per campioni in plastica (polipropilene, PP) da 120 mL, adatto ai rack dei campioni 6.2041.470, 6.2041.840, 6.2041.850, 6.2068.010, 6.2068.020, 6.2068.030 e 6.9920.141, 250 pezzi

### Accessori opzionali

Order no.	Descrizione
6.1458.040	Titration Head, 3 x SGJ 14  Titration head for 120 mL beakers. With 3 openings SGJ 14 and 4 openings 6.4 mm, 2 dosing tips included. Suitable for rinsing and aspiration equipment.



6.1459.307 Sample beaker, plastic (PP), 120 mL, 1 piece


Sample beaker, plastic (polypropylene, PP), 120 mL, suitable for sample racks: 6.2041.470, 6.2041.840, 6.2041.850, 6.2068.010, 6.2068.020, 6.2068.030, and 6.9920.141, 1 piece



---

6.1459.320 Screw cap for 61459300 / 250 pieces

For safe closure of 120 mL sample beakers



---

6.1909.050 Stirring propeller / 102 mm

Propeller for 802 Stirrer when this is used with sample racks for 120 mL beakers (6.1459.300)



---

6.1909.060 Stirring propeller / 95 mm

Stirring propeller with 3 wings for intensive stirring



6.2037.080

Foil holder

Foil holder for 120 mL PP beaker (6.1459.300)



---

6.2041.470

Sample rack, 22 x 120 mL, (PP)

Sample rack for 810 Sample Changer, 814 USB Sample Processor, 855 Robotic Titrosampler "Light", suitable for 22 sample beakers. The following sample beakers can be used: 6.1459.300, 6.1459.307.

Plastic: Polypropylene (PP)



---

6.2041.840

Sample rack, 59 x 120 mL, (PVC)

Sample rack for 815 Robotic USB Sample Processor XL, 855 Robotic Titrosampler "Basic", suitable for 59 sample beakers. The following sample beakers can be used: 6.1459.300, 6.1459.307.

Plastic: Polyvinyl chloride (PVC)



---

6.2041.850

Sample rack, 59 x 120 mL, (PP)

Sample rack for 815 Robotic USB Sample Processor XL, 855 Robotic Titrosampler "Basic", suitable for 59 sample beakers. The following sample beakers can be used: 6.1459.300, 6.1459.307.


Plastic: Polyvinyl chloride (PVC)



6.2068.010      Sample rack, 20 x 120 mL, (PP)

Sample rack for 864 Balance Sample Processor, suitable for 20 sample beakers. The following sample beakers can be used: 6.1459.300.

Plastic: Polypropylene (PP)



---

6.2068.020      Sample rack, 24 x 120 mL, (PVC)

Sample rack for 815 Robotic Soliprep, suitable for 24 sample beakers. The following sample beakers can be used: 6.1459.300, 6.1459.307.

Plastic: Polyvinyl chloride (PVC)



---

6.9920.141      Sample rack, 53 x 120 mL, (PP)

Sample rack for 815 Robotic USB Sample Processor XL, 855 Robotic Titrosampler "Basic", suitable for 53 sample beakers. The following sample beakers can be used: 6.1459.300, 6.1459.307.

Plastic: Polypropylene (PP)

