



Rack dei campioni, 24 x 75 mL, (PVC)

6.2041.340

Rack dei campioni per 810 Sample Processor, 814 USB Sample Processor, 855 Robotic Titrosampler "Light" adatto a 24 beaker per campioni. Possono essere utilizzati i seguenti beaker per campioni: 6.1432.210, 6.1459.400.

Plastica: cloruro di polivinile (PVC)

Accessori opzionali

Order no.	Descrizione
6.1432.210	Sample beaker, clear glass, 75 mL, 1 piece
	Sample beaker, clear glass, 75 mL, suitable for sample racks: 6.2041.340, 6.2041.350, 6.2041.800 and 6.2041.900, 1 piece



6.1458.020

Titration head / 4 x M10 connectors

With 4 M10 screw nipples, including stoppers.



6.1458.030

KF titration head / 4 x M10 connectors

With 4 M10 screw nipples, including stoppers.



6.1458.070

Titration head / 2 x SGJ 14

Titration head for 75 mL beakers. With 2 openings SGJ 14, 2 openings 6.4 mm, and 3 pre-mounted dosing tips. Suitable for rinsing and aspiration equipment.



6.1459.400

Sample beaker, plastic (PE), 75 mL, 144 pieces

Sample beaker, plastic (polyethylene, PE) 75 mL, suitable for sample racks: 6.2041.340, 6.2041.350, and 6.2041.800, suitable for 6.02041.040 OMNIS Software sample rack, 144 pieces.

