



TX-900MH

2.143.0200

USBRS232

- 900 Touch Control (+)
- 915 KF Ti-Touch (+)
- 916 Ti-Touch (+)
- 917 Coulometer (+)
- 91X Meter (+)
- Eco Dosimat / Titrator (+)
- 877 / 848 Titrino plus (6.2151.100)
- 865 / 876 Dosimat plus (6.2151.100)
- 862 Compact Titrosampler (6.2151.100)
- 870 KF Titrino plus (6.2151.100)
- 899 Coulometer (6.2151.100)

(+)

2.143.0200

Qt.	Order no.	
1 PCS	1.143.0200	TX-900MH impact printer
Compact impact printer with USB and RS232 interface for 900 Touch Control 915 KF Ti-Touch 916 Ti-Touch 917 Coulometer 877/848 Titrino plus 865/876 Dosimat plus 91X meter (cable 6.2151.140) Eco Dosimat Eco Titrator Eco Coulometer 862 Compact Titrosampler 870 KF Titrino plus 899 Coulometer		



1 PCS	6.2134.040	RS-232 DB9 - IBM-PC(DB9) connecting cable
-------	-------------------	--

For the Metrohm RS-232 interface - PC connection. For Titrinos, Titrino plus and 756, 831 KF Coulometer, RS-232/USB Box



1 PCS	6.2151.020	Cable USB A - USB B 1.8 m
-------	-------------------	----------------------------------

USB connecting cable



Y cable for connection of a USB printer to the pH/Conductometers 912 / 913 / 914. This cable allows a printer and the power supply unit to be connected to the measuring instrument at the same time.



Order no.

6.2151.100 Adapter USB MINI (OTG) - USB A

For connecting USB instruments.



6.9001.070 Paper roll for TX-900MH

Paper roll for TX-900MH impact printer



6.9001.080 Continuous adhesive labels for TX-900MH

Self-adhesive paper roll for TX-900MH impact printer



6.9002.070 Color ribbon for TX-900MH impact printer

