



804 Ti Stand without stand rod

2.804.0010

Titration stand and controller for 802 Stirrer. The 804 Ti Stand together with the optional 802 Stirrer provides an alternative to the magnetic stirrer. Without stand.

2.804.0010

Qt.	Order no.	
1 PCS	1.804.0010	804 Ti Stand without stand rod
	Titration stand and controller for 802 Stirrer. The 804 Ti Stand together with the optional 802 Stirrer provides an alternative to the magnetic stirrer. Without stand.	




Order no.

6.2001.090

Base plate on both sides

Double stand for attaching 2 x 801, 803, 804 to a Titrande



6.2151.010

Cable MSB plug / socket

Extension cable for MSB connections and for Touch Control - Titrande,
Mini-DIN socket - mini-DIN plug

