



# Pt1000 F

6.00500.600

Pt1000

:

- $\text{F}^-$  ( $10^{-6}$  mol/L)
- 
- 

():  $c(\text{KCl}) = 3 \text{ mol/L}$

()

Order no.

**6.02104.600      Electrode cable plug-in head U / plug P, 0.55 m**

Cable for connecting an electrode with plug-in head U to an analog measuring module from OMNIS.



---

**6.02104.610      Electrode cable plug-in head U / plug P, 1.5 m**

Cable for connecting an electrode with plug-in head U to an analog measuring module from OMNIS.



---

**6.2104.600      Electrode cable for plug in head U/plug F, 2x2 mm B, 1m**

For connecting electrodes with Metrohm plug-in head U to Metrohm instruments (socket F).



---

**6.2104.610      Electrode cable for plug in head U/plug F, 2x2 mm B, 2m**

For connecting electrodes with Metrohm plug-in head U to Metrohm instruments (socket F).



6.2308.020

**Electrolyte 3 mol/L KCl (250 mL)**

Electrolyte solution  $c(\text{KCl}) = 3 \text{ mol/L}$ , (for Ag/AgCl reference systems)

