



Metrosep A Supp 17 - 150/4.0

6.01032.420

Metrosep A Supp 17 - 150/4.0 1.4 mL/min Metrosep A-Supp-17

6.01032.420

Qt.	Order no.	
2 PCS	6.2744.060	Threaded stopper
For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.		



Order no.**6.01032.500 Metrosep A Supp 17 Guard/4.0**

The Metrosep A Supp 17 Guard/4.0 reliably protects the Metrosep A Supp 17 analytical separation columns against contamination. Thanks to the "On Column Guard System", the guard column is very easy to handle. The guard column screws easily onto the analytical column. No tools are required.



6.01032.510 Metrosep A Supp 17 S-Guard/4.0

The Metrosep A Supp 17 S-Guard/4.0 reliably protects the Metrosep A Supp 17 analytical separation columns against contamination. The guard column is simply connected to the separation column using capillary connections. No tools are required.



6.01032.530 Metrosep A Supp 17 S-Guard - 50/4.0

The Metrosep A Supp 17 S-Guard 50/4.0 reliably protects the Metrosep A Supp 17 analytical separation columns against contamination. The guard column is simply connected to the separation column using capillary connections. No tools are required.

