



dThermoprobe HF

6.01118.300

HFOMNISTET
Thermoprobe
:

-
-

6.01118.300

Qt.	Order no.	
1 PCS	6.1236.050	Sleeve with SGJ 14/12 mm
	Sleeve with SGJ 14/12 mm, polyethylene.	



Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Order no.

6.02104.300 Electrode cable plug-in head Q / plug P, 0.55 m

Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.



6.02104.310 Electrode cable plug-in head Q / plug P, 1.5 m

Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.

