



## 75 mL30

6.01402.003

75 mL OMNIS 602041040 30

### Order no.

**6.01403.010      Titration head 4xM10 (P&P)**

Titration head with 4 M10 openings for use in Pick&Place modules of the OMNIS Sample Robot Pick&Place



**6.01403.020 Titration head 4xM10 for KFT (P&P)**

Titration head with 4 M10 openings for Karl Fischer titration in Pick&Place modules of the OMNIS Sample Robot Pick&Place with 4 M10 screw nipples and stoppers.



---

**6.01403.040 Titration head 2xNS14 (P&P)**


Titration head with 2 NS14 openings for use in Pick&Place modules of the OMNIS Sample Robot Pick&Place with 3 retracted dosing tips. Prepared for being equipped with rinsing nozzles. One aspiration tube is already retracted into the titration head.



---

**6.01404.040 Beaker adapter for OMNIS 75 mL sample beaker**

Beaker adapter for centering the 75 mL sample beaker in the slide of the Pick&Place modules



---

**6.01900.020 Stirring propeller 13 mm ETFE**

Stirring propeller 13 mm ETFE with adjustable height for use with Rod Stirrer "Titrator" or Rod Stirrer "Sample Robot"




**6.02041.040**

**OMNIS sample rack, 25 x 75 mL, (PP)**

OMNIS sample rack for OMNIS Sample Robot Pick&Place, suitable for 25 sample beakers. The following sample beakers can be used: 6.01402.000, 6.01402.003, 6.1459.400.

Plastic: Polypropylene (PP)



---

**6.02710.000**

**KF Dis-Cover lid for OMNIS 75 mL sample beaker, 25 pieces**

Lid with seals for volumetric Karl Fischer Titration in 6.01402.000 and 6.01402.003 sample beakers in the OMNIS Sample Robot Pick&Place system.



---

**6.02710.030**

**Dis-Cover lid for OMNIS 75 mL sample beaker, 25 pieces**

Dis-Cover lid for 75 mL 6.01402.000 and 6.01402.003 sample beakers in the OMNIS Sample Robot Pick&Place system.

