



3xNS14/ 4x6.4mm (P&P)

6.01403.030

2OMNIS Sample Robot Pick&Place Pick&Place 120 mL 3 NS14 4 x 6.4 mm

6.01403.030

Qt.	Order no.	
2 PCS	6.1236.020	Sleeve with SGJ 14/12 mm
	Sleeve with SGJ 14/12 mm and O-ring.	



2 PCS

6.1446.000

SGJ stopper / B-14/(15)



2 PCS

6.1808.000

Connection bushing M6



2 PCS

6.2726.090

Antidiffusion valve

For antidiffusion buret tips 6.1541.XXX, 6.1543.130, 6.1543.140,
6.1543.150, 6.1543.160, 6.1543.200



Order no.

6.1829.050

Replacement tubings for retracted tubings in titration heads and robotic arms

Three replacement tubings with tapering tip.

