



13 mm ETFE

6.01900.020

Titration Sample Robot ETFE 13 mm

6.01900.020

Qt.	Order no.	
0 PCS	6.2803.000	Grease (1 oz.)
	Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX	



Order no.

6.01402.000 Sample beakers, clear glass, 75 mL, 30 pieces

Sample beakers, clear glass, 75 mL, suitable for 602041040 OMNIS sample rack, 30 pieces



6.01403.040 Titration head 2xNS14 (P&P)

Titration head with 2 NS14 openings for use in Pick&Place modules of the OMNIS Sample Robot Pick&Place with 3 retracted dosing tips. Prepared for being equipped with rinsing nozzles. One aspiration tube is already retracted into the titration head.

