



OMNIS9 x 250 mL(PP)

6.02041.010

9 OMNIS& OMNIS: 6.01400.0006.01400.0036.01400.100

: (PP)

6.02041.010

Qt.	Order no.	
1 PCS	6.02042.000	Holder for rack label (3x)

The holder for rack labels permits durable labeling of the rack with an unambiguous name.



Order no.

6.01400.000 Sample beakers, clear glass, 250 mL, 10 pieces

Sample beakers, clear glass, 250 mL, suitable for 602041010 OMNIS sample rack, 10 pieces.



6.01400.003 Sample beakers, amber glass, 250 mL, 10 pieces

Sample beakers, amber glass, 250 mL, suitable for 602041010 OMNIS sample rack, 10 pieces.



6.01400.100 Sample beakers, plastic (PP), 250 mL, 10 pieces

Sample beakers, plastic (polypropylene, PP), 250 mL, suitable for 602041010 OMNIS sample rack, 10 pieces.

