



6.02743.090

PP 50 mL (6.02743.030) PE 125 ml (6.1627.000) PE 50 ml (6.1608.100)

858 Professional Sample Processor (6.2624.200) Pick-up (6.2624.100)

**Order no.**

**6.02743.030      Bottle, PP, 50 mL, 12 pieces**

The 50 mL bottles are made of polypropylene (PP) plastic and can be used for samples, standards, check standards and rinsing solutions. Production and packaging both take place under sterile and clean room conditions. This guarantees minimum contamination. The sample vessels are subjected to random checks for leachable anions and cations.



**6.1608.100**

**PE bottle 50 mL**

Bottle for Extension Module Liquid Handling.



---

**6.1627.000**

**PE bottle, 125 mL for IC**

This bottle is used with fully automatic dilution.



---

**6.2624.100**

**Pick-up needle MF**

Coated steel needle for metal-free pick-up sampling as well as for Eppendorf sample vials.



---

**6.2624.200**

**Sample needle made of steel, coated**

Alternative to the sample needle made of zirconium oxide. Particularly suitable for septum vessels. Metal-free thanks to special coating. Requires the 1/16 in. needle holder (6.2833.010) for installation on the Sample Processor.

