



6.05800.030

OMNIS Sample Robot Pick&Place (2 ) (4 )

## 6.05800.030

Qt.	Order no.	
1 PCS	6.01543.100	Aspiration tip automation

Aspiration tip for automated OMNIS titrations.




3 PCS

**6.01805.000**

**FEP tubing connection / M6 / 66.5 cm**

With light protection and protection against kinks.



---

1 PCS

**6.01805.010**

**PTFE tubing connection / M8 / 66.5 cm**

With light protection and protection against kinks.



---

1 PCS

**6.05001.010**

**"Peristaltic" (2-channel) pump module  
accessory set**

The accessory set contains all of the components required for the tube connections between a rinsing pump and an aspiration pump at the "Peristaltic" pump module. In addition to preassembled tubing with Luer connectors, each of these also includes two press clamps and pump tubing.



---

2 PCS

**6.1621.000**

**Container 10 L**

As rinsing or waste container in automated systems.



2 PCS

**6.1621.100**

**Cap for 10-L container**

Cap for 10-L container. For easy and safe transportation of filled containers.



---

2 PCS

**6.1828.000**

**PVDF connection nipple**

For 6.1621.000 canister.



---

1 PCS

**6.2709.070**

**Guiding sleeve**

For fixation of buret and dosing tips in titration heads and electrode holders



---

3 PCS

**6.2740.020**

**Spray nozzle**

For the fine-spraying of the rinsing solution.

