



Metrosep C 3 - 250/4.0

6.1010.430

Metrosep C 3 - 250/4.042'00051'00031'000

Metrosep C 3 - 250/4.0Metrosep C 3 - 250/4.0

Order no.

6.1010.450

Metrosep C 3 Guard/4.0

The Metrosep C 3 S-Guard/4.0 contains the C-3 column material and is used to protect Metrosep C-3 cation columns. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column. The Metrosep C 3 Guard/4.0 also works based on the "On Column Guard System" and is screwed directly onto the separation column with virtually no dead volume.



The Metrosep C 3 S-Guard/4.0 contains the C-3 column material and is used to protect Metrosep C-3 cation columns. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column.

