



Hamilton RCX-30 - 250/4.6

6.1018.000

Hamilton RCX-30 - 250/4.6 -RCX-30 - 250/4.6

Hamilton RCX-30 - 250/4.6 2 mL/min (250 mm) ()

Order no.

6.1011.030 Metrosep RP 2 Guard/3.5

The Metrosep RP 2 Guard/3.5 is a guard column for universal use. It reliably protects the analytical separation column against contaminations. It reliably protects the analytical separation column against contamination, removing the smallest particles, traces of iron oxide and bacteria. The Metrosep RP 2 Guard/3.5 helps to reduce costs; its filter disk can be replaced very easily.



6.1011.040

Metrosep RP 3 Guard HC/4.0

The Metrosep RP 3 Guard HC/4.0 is a guard column for universal use. It reliably protects the analytical separation column against contamination, securely removing lipophilic organic contamination, the smallest particles, traces of iron oxide and bacteria. The guard column is based on a polymer material and thanks to its larger pack volume, has a significantly higher capacity than the Metrosep RP 2Guard/3.5. It can be used throughout the entire pH range.



6.1011.130

Spare filter for RP 2 Guard/3.5, 10 pieces

