



Metrosep A Supp 16 - 250/2.0

6.1031.230

Metrosep A Supp 16 - 250/2.0 () MSM-LCMetrosep A Supp 162mm4mm/Metrosep A Supp 16 - 250/2.0

Metrosep A Supp 16 - 250/2.0 IC-MS

6.1031.230

Qt.	Order no.	
2 PCS	6.2744.060	Threaded stopper
For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.		



Order no.**6.1031.600****Metrosep A Supp 16 Guard/2.0**

The Metrosep A Supp 16 Guard/2.0 reliably protects the Metrosep A Supp 16 analytical separation columns with 2 mm inner diameter against contamination. Thanks to the "On Column Guard System", the guard column is very easy to handle. The guard column screws easily onto the analytical column. No tools are required.

**6.1031.610****Metrosep A Supp 16 S-Guard/2.0**

The Metrosep A Supp 16 S-Guard/2.0 reliably protects the Metrosep A Supp 16 analytical separation columns with 2 mm inner diameter against contamination. The guard column is connected to the separation column simply using capillary connections. No tools are required.

