



Buffer solution pH 7.00 (30 x 30 mL)

6.2307.210

Buffer solution pH 7.00 (25 °C) in single use sachets, colorless

Order no.

6.2307.100 Buffer solution pH 4 (500 mL)

Ready-to-use buffer solution pH = 4 (25 °C), colored, with closure seal



6.2307.110

Buffer solution pH 7 (500 mL)

Ready-to-use buffer solution pH = 7 (25 °C), colored, with closure seal



6.2307.120

Buffer solution pH 9 (500 mL)

Ready-to-use buffer solution pH = 9 (25 °C), colorless, with closure seal



6.2307.200

Buffer solution pH 4.00 (30 x 30 mL)

Buffer solution pH 4.00 (25 °C) in single use sachets, colorless



6.2307.220

Buffer solution pH 9.00 (30 x 30 mL)

Buffer solution pH 9.00 (25 °C) in single use sachets, colorless



6.2307.230

Buffer solutions pH 4, 7 and 9

Mixed buffer solutions pH 4.00/7.00/9.00 (25 °C) in single use sachets, colourless, box of 3 x 10 x 30 mL



6.2325.000

pHit kit

Maintenance kit for pH electrodes

The kit contains:

- 50 mL cleaning solution
- 50 mL 3M KCl solution
- 50 mL storage solution
- 2 Storage vessels
- Instructions for use

