



6.2621.080

1/16" 1/8" 5

**Order no.**

**6.1803.040**      **PTFE capillary 0.5 mm i.d. / 1 m**

Capillary for sample handling in IC.



**6.1803.080**

**PTFE capillary 0.97 mm i.d. / 1 m**

Capillary for filtration, for IC Filtration Sample Processor, IC Liquid Handling Filtration Unit



---

**6.1831.010**

**PEEK capillary 0.25 mm i.d. / 3 m**

For all IC components.



---

**6.1831.030**

**PEEK capillary 0.75 mm i.d. / 3 m**

For all IC components.



---

**6.1831.060**

**PEEK capillary 0.5 mm i.d., 1 m**

For inline dialysis and for Sample Changers. For IC Dialysis Unit, IC Sample Processor, Compact IC Autosampler, IC Filtration Sample Processor, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit

