



Low Volume

6.2729.200

4 6.2744.060

6.2729.200

Qt.	Order no.	
4 PCS	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Order no.

6.2714.010 Dialysis membrane (cellulose acetate)

Nominal pore size : 0.2 µm. Set of 50 pieces. For inline dialysis



6.2714.030 Dialysis membrane (polyamide)

For inline dialysis and filtration



6.2744.000 Pressure screw PVDF

With UNF 10/32 connection. For IC instruments, inline dialysis (connection of PTFE capillaries)

