



## Accessory kit: Vario/Flex ChS

6.5000.030

Set of accessories for chemical suppression (ChS) for Professional IC Vario and Compact IC Flex instruments. Without Suppressor Rotor.

**6.5000.030**

Qt.	Order no.	
1 PCS	<b>6.1602.150</b>	<b>Bottle attachment / GL 45 - 3 x UNF 10/32</b>

For the connection of 1/16 in. capillaries. Used in conjunction with MSM auxiliary solutions and in Metrohm Inline Sample Preparation.



1 PCS

**6.1608.020**

**Glass bottle / 1000 mL / GL 45**

Bottle for auxiliary solutions.



1 PCS

**6.1803.030**

**PTFE capillary 0.5 mm i.d. / 3 m**

Capillary for inline dialysis, for Dialysis Unit, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit



2 PCS

**6.2744.040**

**2 x UNF 10/32 coupling**

For connecting 1/16 in. capillaries. For IC instruments



1 PCS

**6.2744.070**

**Pressure screw short**

Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries



Universal connecting piece for MSM in Professional IC Vario and Compact Flex instruments and for MSM-HC and SPM module. With inlets and outlets.

