



826 pH mobile

2.826.0010

Portable pH meter with IrDa interface, including carrying strap and batteries.

Phạm vi giao hàng 2.826.0010

Qt.	Order no.	Description
1 PCS	1.826.0010	826 pH mobile

Portable pH meter with IrDa interface, including carrying strap and batteries.



1 PCS

6.2050.000

Carrying strap for pH Meter

Carrying strap for 704, 826 pH meter.



1 PCS

6.2133.000

Batteries LR6, 1.5V, 4 pieces

1.5 V LR6 type. Set of 4 pieces. Used with 704 pH Meter, 773 pH/mV Simulator



1 PCS

6.2621.140

Hexagon key 2.5 mm



Phụ kiện tùy chọn

Order no.	Description	
6.2307.230	Buffer solutions pH 4, 7 and 9 Mixed buffer solutions pH 4.00/7.00/9.00 (25 °C) in single use sachets, colourless, box of 3 x 10 x 30 mL	
6.2325.000	pHit kit Maintenance kit for pH electrodes The kit contains: • 50 mL cleaning solution • 50 mL 3M KCl solution • 50 mL storage solution • 2 Storage vessels • Instructions for use	