



20 – 90 mL titration vessel

6.01406.220

20 – 90 mL titration vessel with screw top for use with 6014050X0 titration vessel lids.

Phụ kiện tùy chọn

Order no.	Description
6.01405.000	5 x SGJ 14 titration vessel lid OMNIS
	Titration vessel lid with five SGJ 14 openings for 601406XX0 titration vessels with screw cap.



6.01405.010 OMNIS titration vessel lid for volumetric Karl Fischer titration

Titration vessel lid for volumetric Karl Fischer titration with an OMNIS titrator.



6.01405.040 OMNIS titration vessel lid for 601406XXX titration vessels for Polytron

Titration vessel lid for volumetric Karl Fischer titration with an OMNIS titrator with Polytron.

