



Bottle cap, multi-use

6.01601.000

Reusable bottle cap with 3S technology for connecting reagent bottles to an OMNIS Titrator, Titration Module or Dosing Module.

Phạm vi giao hàng 6.01601.000

Qt.	Order no.	Description
1 PCS	6.02701.010	PTFE sealing ring for 601601000 (5 pcs)

PTFE sealing ring for bottle cap, multi-use.



1 PCS

6.1819.020

Tube

For M6 threaded holes. For 6.1602.1XX bottle attachments

