



LL-Ag/AgCl reference electrode for CVS

6.0728.130

Inner reference system with electrolyte $c(\text{KCl}) = 3 \text{ mol/L}$. For use with electrolyte vessel made of plastic 6.1245.010.

Phụ kiện tùy chọn

Order no.	Description
6.1245.010	Electrolyte vessel made of plastic
	Electrolyte vessel made of plastic with ceramic diaphragm for double junction reference electrode.



6.2112.130

Adapter plug contact 5.5 mm / screw contact, gold-plated

For VA electrodes

